



D-Amphetamine [KLH] (DAGA-244K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	D-Amphetamine (Dextroamphetamine) is a central nervous system stimulant and a sympathomimetic. It has also been used in the treatment of narcolepsy and of attention deficit disorders and hyperactivity in children.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords D-Amphetamine; Dextroamphetamine