



Cyclopentobarbital [HRP] (DAGA-246H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Barbiturates are a class of around 12 compounds derivatised from barbituric acid.
Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Cyclopentobarbital
----------	--------------------