



Dapsone [OVA] (DAGA-248O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Dapsone , a sulfone, belongs to the family of medicines called anti-infectives. Dapsone is used to treat leprosy (Hansen's disease) and to help control dermatitis herpetiformis, a skin problem. When it is used to treat leprosy, dapsone may be given with one or more other medicines.
--------------	--