



Dehydronorketamine [KLH] (DAGA-250K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Dehydronorketamine is a major metabolite of Ketamine which is an NMDA receptor antagonist has a wide range of effects in humans, including analgesia, anesthesia, hallucinations, elevated bp, and bronchodilation.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Dehydronorketamine
-----------------	--------------------