



Digoxigenin [BSA] (DAGA-257B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Digoxigenin (DIG) is a steroid found exclusively in the flowers and leaves of the plants <i>Digitalis purpurea</i> and <i>Digitalis lanata</i> . Digoxigenin is chemically closely related to Digoxin, the cardiac glycoside used for the treatment of various heart diseases.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	DIG; Lanadigenin;
-----------------	-------------------