



## Digoxigenin [KLH] (DAGA-257K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Digoxigenin (DIG) is a steroid found exclusively in the flowers and leaves of the plants <i>Digitalis purpurea</i> and <i>Digitalis lanata</i> . Digoxigenin is chemically closely related to Digoxin, the cardiac glycoside used for the treatment of various heart diseases.
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** DIG; Lanadigenin;