



Dipyron [BSA] (DAGA-258B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Dipyron or "Mexican aspirin" is a nonsteroidal anti-inflammatory drug. The mechanism of action of dipyron is thought to be similar to that of other NSAIDs: inhibition of the production of prostaglandins.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Analgin; C13H16N3O4NaS; Conmel; Metamizole sodium; Mexican aspirin; Neo Melubrina; Prodolina F; Prodolina F;
-----------------	--