



ONNV Capsid protein [His] (DAGA-260)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The o'nyong'nyong virus or o'nyong-nyong virus (ONNV) is a virus first isolated by the Uganda Virus Research Institute in Entebbe, Uganda in 1959. It is a togavirus (family Togaviridae), genus Alphavirus and is closely related to the chikungunya and Igbo
Purity	>90% , based on SDS PAGE
Conjugate	His
Applications	WB standard, antibody ELISA, immunogen, etc.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg, 1 mg
Buffer	In PBS with 8M Urea
Preservative	None
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.

BACKGROUND

Keywords	o'nyong-nyong virus capsid protein;ONNV capsid protein;ONNV C protein;o'nyong'nyong virus protein C;ONNV
-----------------	--