



Eddp [BSA] (DAGA-260B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords EDDP; EDIFEN; EDIFENFOS; EDIFENPHOS; EDISAN; edinphensoph; EDIFENPHOS PESTANAL (O-ETHYL S; S-DIPHENY); EDIFENPHOS; 500MG; NEAT; eddp (jmaf); edifenphos (bsi; iso); ENDINFENPHOS; O-ETHYLDIPHENYLPHOSPHORODITHIOLATE; EDINFENPHOS; edifenphos (ISO) O-ethyl S; S-diphenyl phosphorodithioate;