



RRV Envelope Antigen 2 (aa 271-334) [His] (DAGA-262)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Ross River virus belongs to the virus genus Alphavirus, which is part of the family Togaviridae. The alphaviruses are a group of small enveloped single-strand positive-sense RNA viruses.
Purity	>95% , based on SDS PAGE
Conjugate	His
Applications	WB standard, antibody ELISA, immunogen, etc.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg, 1 mg
Buffer	In PBS with 8M Urea
Preservative	None
Storage	Stable for 6 months at -20°C from the date of shipment. Nonhazardous. No MSDS required.

BACKGROUND

Keywords	RRV E2;Ross River Virus envelope protein;RRV Envelope Antigen 2;RRV
----------	---