



Escitalopram [KLH] (DAGA-262K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Escitalopram is an antidepressant of the selective serotonin reuptake inhibitor (SSRI) class. It is used in the treatment of adults with major depressive disorder, generalized anxiety disorder, social anxiety disorder, or panic disorder.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	(S)-(+)-CITALOPRAM OXALATE; (1R)-1-(3-dimethylaminopropyl)-1-(4-fluorophenyl)-3H-isobenzofuran-5-carbonitrile; 1-(3-Dimethylaminopropyl)-1-(4-fluorophenyl)-1, 3-dihydroisobenzofuran-5-carbonitrile; Escitalopram; Anxiset E(India) Lexapro; Cipralex; Seroplex; Lexamil; Lexam; Lexapro
-----------------	---