



Ethylglucuronide [KLH] (DAGA-263K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Ethylglucuronide is a metabolite which is formed in the body following ingestion of ethanol. It is used as a biomarker to test for ethanol use in situations where such use is prohibited, such as by the military or in recovering alcoholic patients.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	EtG;
-----------------	------