



Halofuginone [KLH] (DAGA-270K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Halofuginone is a coccidiostat used in veterinary medicine. It is a synthetic halogenated derivative of febrifugine, a natural quinazolinone alkaloid which can be found in the Chinese herb Dichroa febrifuga. Halofuginone inhibits the development of T helper28cells, immune cells that play an important players in autoimmune disease, but it does not affect other kinds of T cells which involved in normal immune function. Halofuginone therefore has potential for the treatment of autoimmune disorders. Halofuginine is also an inhibitor of collagen type I gene expression and as a consequence it may inhibit tumor cell growth.
----------------------------	--

Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Halofuginone; HF
-----------------	------------------