



## Haloperidol [KLH] (DAGA-271K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Haloperidol is a butyrophenone antipsychotic drug that has strong central antidopaminergic activity. It is used in the treatment of schizophrenia, acute psychotic states and delirium.
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

<b>Keywords</b>	Aloperidin; Bioperidolo; Brotopon; Butyrophenone antipsychotic drug; Dozic; Einalon S; Eukystol; Haldol; Halosten; Keselan; Linton; Peluces; Serenace; Serenase; Sigaperidol;
-----------------	---