



## Human PCP4 [HRP] (DAGA-272H)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                                                                                       |
|---------------|---------------------------------------------------------------------------------------|
| Conjugate     | HRP                                                                                   |
| Applications  | ELISA                                                                                 |
| Format        | Liquid                                                                                |
| Concentration | Batch dependent - please inquire should you have specific requirements.               |
| Size          | 1 mg                                                                                  |
| Buffer        | 0.01M pH7.4 PBS                                                                       |
| Preservative  | None                                                                                  |
| Storage       | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

### BACKGROUND

|          |                                                                                                                                                   |
|----------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Keywords | PCP4; Purkinje cell protein 4; PEP 19; brain-specific antigen PCP-4; brain specific polypeptide PEP19; brain-specific polypeptide PEP-19; PEP-19; |
|----------|---------------------------------------------------------------------------------------------------------------------------------------------------|