



## Human PCP4 [KLH] (DAGA-272K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

Keywords	PCP4; Purkinje cell protein 4; PEP 19; brain-specific antigen PCP-4; brain specific polypeptide PEP19; brain-specific polypeptide PEP-19; PEP-19;
----------	---