



Imipramine [KLH] (DAGA-275K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Imipramine (sold as Antidepressin, Deprimin, Deprinol, Depsol, Depsonil, Dynaprin, Eupramin, Imipramil, Irmin, Janimine, Melipramin, Surplix, Tofranil), also known as melipramine, is an antidepressant medication, a tricyclic antidepressant of the dibenzazepine group. Imipramine is mainly used in the treatment of major depression and enuresis (inability to control urination).
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Antidepressin; Deprimin; Deprinol; Depsol; Depsonil; Dynaprin; Eupramin; Imipramil; Irmin; Janimine; Melipramin; Surplix; Tofranil; melipramine; Imipramine
-----------------	---