



## JWH 018 [HRP] (DAGA-277H)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Antigen Description	JWH 018 is a mildly selective agonist of the peripheral cannabinoid (CB2) receptor, derived from the aminoalkylindole WIN 55,212-2.
Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** JWH 018; 1-pentyl-3-(1-naphthoyl)indole); AM678;