



JWH-250 [HRP] (DAGA-278H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords 1-pentyl-3-(2-methoxyphenylacetyl)indole; 2-(2-methoxyphenyl)-1-(1-pentylindol-3-yl)ethanone;