



## Ketoprofen [BSA] (DAGA-281B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Ketoprofen, (RS)2-(3-benzoylphenyl)-propionic acid, is a common non-steroidal anti-inflammatory drug (NSAID) with analgesic and antipyretic effects. It acts by inhibiting the body's production of prostaglandin
<b>Conjugate</b>	BSA
<b>Applications</b>	LFIA
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** Ketoprofen; 2-(3-benzoylphenyl)-propionic acid; ORUDIS; ORUVAIL;