



## Ketoprofen [KLH] (DAGA-281K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Antigen Description	Ketoprofen, (RS)2-(3-benzoylphenyl)-propionic acid, is a common non-steroidal anti-inflammatory drug (NSAID) with analgesic and antipyretic effects. It acts by inhibiting the body's production of prostaglandin
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** Ketoprofen; 2-(3-benzoylphenyl)-propionic acid; ORUDIS;