



## Lasalocid A [KLH] (DAGA-283K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Lasalocid is able to make neutral complexes with monovalent and divalent cations and transport them through apolar phase (including lipid bilayer membranes). It can transport also big organic cations like dopamine.
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

<b>Keywords</b>	Lasalocid A
-----------------	-------------