



Levomethadyl acetate [BSA] (DAGA-285B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Levomethadyl acetate is indicated as a second-line treatment for the treatment and management of opioid dependence if patients fail to respond to drugs like methadone or buprenorphine.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Levomethadyl acetate; Levo Alpha Acetylmethadol; Orlaam;
----------	--