



Lidocaine [KLH] (DAGA-286K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Lidocaine belongs to the family of medicines called local anesthetics. When lidocaine is applied to the skin, it produces pain relief by blocking the signals at the nerve endings in the skin.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	2 Diethylamino N 27dimethylphenyl acetamide; Acetamide 2-(diethylamino)-N-(26-dimethylphenyl);
----------	--