



Recombinant Human Insulin (DAGA-289)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	≥98% as determined by HPLC
Conjugate	Unconjugated
Cas.No	11061-68-0
Applications	Cell culture, Animal modeling
Molecular Weight	5.8 kDa
Reconstitution	Soluble in water at ≥5 mg/mL (ultrasonication and the addition of HCl may be required to aid dissolution)
Bio-activity	≥27IU/mg
Format	Lyophilized
Size	10 mg, 20 mg, 50 mg, 100 mg
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	Insulin; INS
----------	--------------