



MDMA [KLH] (DAGA-291K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	MDMA (3,4-Methylenedioxymethamphetamine)— colloquially known as ecstasy, often abbreviated "E" or "X"— is an entactogenic drug of the phenethylamine and amphetamine class of drugs. MDMA can induce euphoria, a sense of intimacy with others, and diminished anxiety. Many studies, particularly in the fields of psychology and cognitive therapy, have suggested that MDMA has therapeutic benefits and facilitates therapy sessions in certain individuals, a practice for which it had formally been used in the past.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords MDMA; methylenedioxymethamphetamine;