



## Mescaline [KLH] (DAGA-295K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Mescaline is a naturally-occurring psychedelic alkaloid. It is used primarily as a recreational drug and is also used to supplement various types of meditation and psychedelic therapy. Users typically experience visual hallucinations and radically altered states of consciousness, often experienced as pleasurable and illuminating but occasionally are accompanied by feelings of anxiety or revulsion.
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** Mescaline; 3 4 5-trimethoxyphenethylamine;