



Methadone P [KLH] (DAGA-297K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Antigen Description | Methadone is a synthetic addictive narcotic drug C ₂₁ H ₂₇ NO used especially in the form of its hydrochloride for the relief of pain and as a substitute narcotic in the treatment of heroin addiction. |
| Conjugate | KLH |
| Applications | Immunogen |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01M pH7.4 PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

Keywords Methadone P; Dolophine;