



Methcathinone [KLH] (DAGA-301K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Methcathinone, known on the streets as "Cat," is a structural analogue of methamphetamine and cathinone. Clandestinely manufactured, methcathinone is almost exclusively sold in the stable and highly water soluble hydrochloride salt form.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Methcathinone; 2-(methylamino)propiofenone; 2-Methylamino-1-phenylpropan-one; monomethylpropion; 2-methylamino-1-phenyl-propan-1-one; Methcathinone
-----------------	---