



Methcathinone[BSA] (DAGA-302A)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Methcathinone, known on the streets as "Cat," is a structural analogue of methamphetamine and cathinone. Clandestinely manufactured, methcathinone is almost exclusively sold in the stable and highly water soluble hydrochloride salt form.
Species	N/A
Conjugate	BSA
Applications	This antigen can be used in Lateral flow immunoassays. Note: This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.
Format	Liquid
Concentration	Hapten/Protein: 25/1
Size	1 mg
Buffer	Distilled water
Preservative	0.02% sodium azide
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.