



# Methylprednisolone [KLH] (DAGA-302K)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Methylprednisolone is a synthetic corticosteroid drug with anti-inflammatory properties. It can be used in the treatment of arthritis, severe allergies, asthma and some types of cancer. Methylprednisolone can alter gene expression and inhibit the production of proinflammatory cytokines through a mechanism where it binds specific nuclear receptors. It can also decrease the number of circulating lymphocytes and stimulate apoptosis.
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Keywords** methylprednisolone; Medrol; Medrone; Urbason; Medrate; Dopomedrol;