



Recombinant Norovirus GI (aa 225-520, P Domain) [His] (DAGA-3040)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Purified Recombinant Antigen (P domain 225-520 with His-tag at N-terminus)
Purity	>95% by SDS-PAGE
Conjugate	His
Applications	ELISA, Controls, Vaccine Research
Molecular Weight	30 kDa
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline with 25 mM arginine
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.