



Recombinant Norovirus GI VP1 (Full length) [His] (DAGA-3041)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Purified Recombinant Antigen (Full length VP1 peptide with 6xHis-tag at N-terminus)
Purity	>95% by SDS-PAGE
Conjugate	His
Applications	ELISA, Controls, Vaccine Research
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.