



Recombinant EBV p18 (DAGA-3074)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The E.Coli derived recombinant mosaic protein contains the HHV-4 p18 regions, 1-119 amino acids.
Purity	>95% , based on SDS PAGE
Conjugate	TBD
Applications	WB standard, antibody ELISA, immunogen, etc.
Format	Each vial contains 100 µg of lyophilized protein in 50mM Tris-HCl pH 8.0, 60mM NaCl and 50% glycerol. Reconstitute the protein with 100 µl of Millipore water.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg
Storage	Before reconstitution, stable for 1 year at -20 centigrade from the date of shipment. After reconstitution, stable for a month at 4 centigrade. Nonhazardous. No MSDS required.