



Norbuprenorphine [KLH] (DAGA-308K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|---|
| Antigen Description | Clinically, buprenorphine is utilized as a substitution therapy for opioid dependence and as an analgesic. Buprenorphine is a partial agonist of the mu-opioid receptor. These mu binding sites are discretely distributed in the human brain, spinal cord, and other tissue. The clinical effects of mu receptor agonists are sedation, euphoria, respiratory depression, and and analgesia. |
| Conjugate | KLH |
| Applications | Immunogen |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01M pH7.4 PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

| | |
|-----------------|-------------------------|
| Keywords | Norbuprenorphine; BUPM; |
|-----------------|-------------------------|