



Recombinant Dust Mite Der f 2 [His] (DAGA-3090)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% pure (SDS-PAGE).
Conjugate	His
Applications	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	0.5 mg
Buffer	50 mM MES, pH 6 containing 0.1 M Potassium Chloride.
Preservative	None
Storage	Store at -20°C to -80°C. Aliquot and store at -80°C to avoid multiple freeze/thaw cycles.