



Recombinant Bourbon virus envelope glycoprotein [His] (DAGA-3104)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	6xHis tagged envelope glycoprotein (Bourbon virus)(aa. 19-462, GenBank No. AMN92169).
Purity	≥ 95%
Conjugate	His
Applications	WB standard, antibody ELISA, Functional analysis, etc
Format	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg
Storage	Keep it at 4 centigrade if used within a month. For long term storage, split it into small aliquots and keep at -80 centigrade. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.