



Recombinant Oropouche virus NSs [His, GST] (DAGA-3105)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	N-terminal GST- and C-terminal 6xHis tagged NSs (Oropouche virus) protein (a.a. 1-92)(GenBank No. AJT39492)
Purity	> 90% as determined by SDS-PAGE
Conjugate	His, GST
Applications	WB standard, antibody ELISA, etc
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg
Buffer	PBS with 8M Urea
Storage	Keep it at 4 centigrade if used within a month. For long term storage, split it into small aliquots and keep at -80 centigrade. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.