



Recombinant Oropouche virus Gn [His] (DAGA-3108)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	6xHis tagged Gn (Oropouche virus) protein (a.a. 24-198)(GenBank No. AJT39490)
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	WB standard, antibody ELISA, etc
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg
Buffer	PBS with 8M Urea
Storage	Keep it at 4 centigrade if used within a month. For long term storage, split it into small aliquots and keep at -80 centigrade. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.