



## Horse IgG [FITC] (DAGA-3118)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	Horse Immunoglobulin G (FITC)
<b>Purity</b>	Horse IgG (FITC) was purified by delipidation, salt fractionation, ion exchange chromatography followed by dialysis.
<b>Conjugate</b>	FITC
<b>Applications</b>	IF
<b>Format</b>	Lyophilized from 0.01M Na <sub>3</sub> PO <sub>4</sub> , pH 7.4, with 0.14M NaCl, 10 mg/ml BSA, and 0.01% thimerisol. Immunoglobulin and protease free.
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Storage</b>	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.