



## Mouse IgG1 [HRP] (DAGA-3119)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	Mouse Immunoglobulin G 1 (HRP)
<b>Purity</b>	Mouse IgG1 (HRP) was purified by Protein A affinity chromatography.
<b>Conjugate</b>	HRP
<b>Applications</b>	User optimized
<b>Format</b>	Lyophilized from 0.02M K3PO4, pH 7.2, with 0.15M NaCl, 10 mg/ml BSA, and 0.01% gentamicin sulfate. Immunoglobulin and protease free. DO NOT add NaN3.
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Storage</b>	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.