



Norfentanyl [BSA] (DAGA-311B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Norfentanyl is a metabolite of Fentanyl. It is a potent narcotic analgesic, abuse of which leads to habituation or addiction. It is primarily a mu-opioid agonist. Norfentanyl is also used as adjunct to general anaesthetics, and as an anaesthetic for induction and maintenance.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Norfentanyl
-----------------	-------------