



Bovine IgG (Fab) [FITC] (DAGA-3121)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Bovine Immunoglobulin G (Fab) (FITC)
Purity	Bovine IgG (Fab) (FITC) was purified by delipidation, salt fractionation and ion exchange chromatography followed by dialysis.
Conjugate	FITC
Applications	IF
Format	Lyophilized from 0.02M K3PO4, pH 7.2, with 0.15M NaCl, 10 mg/ml BSA, and 0.01% NaN3. Immunoglobulin and Protease free.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.