



Goat IgG (Fab'2) [HRP] (DAGA-3127)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Immunoglobulin G (Fab'2) (HRP)
Purity	Goat IgG (Fab'2) (HRP) was purified by delipidation, salt fractionation and ion exchange chromatography followed by dialysis.
Conjugate	HRP
Applications	User optimized
Format	Lyophilized from 0.02M K3PO4, pH 7.2, with 0.15M NaCl, 10 mg/ml BSA, and 0.01% gentamicin sulfate. Immunoglobulin and protease free. DO NOT add NaN3.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.