



## Human IgA protein [FITC] (DAGA-3128)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	Human Immunoglobulin A protein (FITC)
<b>Purity</b>	Human IgA protein (FITC) was purified by delipidation, salt fractionation, ion exchange chromatography and dialysis.
<b>Conjugate</b>	FITC
<b>Applications</b>	IF
<b>Format</b>	Lyophilized from 0.02M K3PO4, pH 7.2, with 0.15M NaCl, 10 mg/ml BSA and 0.01% NaN3. Immunoglobulin and protease free.
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Storage</b>	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.