



Human IgA protein [HRP] (DAGA-3133)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Human Immunoglobulin A protein (HRP)
Purity	Human IgA protein (HRP) was purified by delipidation, selective precipitation, tandem molecular sieve chromatography and dialysis.
Conjugate	HRP
Applications	User optimized
Format	Lyophilized from 0.02M K ₃ PO ₄ , pH 7.2, with 0.15M NaCl, and 0.01% gentamicin sulfate. DO NOT add NaN ₃ .
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.