



## Goat IgG (Fab) [HRP] (DAGA-3137)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                         |                                                                                                                                                                     |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b> | Goat Immunoglobulin G (Fab) (HRP)                                                                                                                                   |
| <b>Purity</b>           | Goat IgG (Fab) (HRP) was purified by delipidation, salt fractionation and ion exchange chromatography followed by dialysis.                                         |
| <b>Conjugate</b>        | HRP                                                                                                                                                                 |
| <b>Applications</b>     | User optimized                                                                                                                                                      |
| <b>Format</b>           | Lyophilized from 0.02M K3PO4, pH 7.2, with 0.15M NaCl, 10 mg/ml BSA, and 0.01% gentamicin sulfate. Immunoglobulin and protease free. DO NOT add NaN3.               |
| <b>Concentration</b>    | Batch dependent - please inquire should you have specific requirements.                                                                                             |
| <b>Size</b>             | 1 mg                                                                                                                                                                |
| <b>Storage</b>          | Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles. |