



Human IgG Fc Fragment (FITC) (DAGA-3138)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Human Immunoglobulin G Fc
Purity	Human IgG Fc was purified by ion exchange chromatography followed by papain digestion and extensive dialysis against the buffer
Conjugate	FITC
Applications	IF, WB
Format	Supplied in lyophilized form in 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2. Restore with deionized water
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 2-8 centigrade.