



Mouse IgG2b [HRP] (DAGA-3139)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse IgG2b (HRP)
Purity	Mouse IgG2b (HRP) was purified by Protein A affinity chromatography.
Conjugate	HRP
Applications	User optimized
Format	Lyophilized from 0.02M K3PO4, pH 7.2, with 0.15M NaCl, 10 mg/ml BSA, and 0.01% gentamicin sulfate. Immunoglobulin and protease free. DO NOT add NaN3.
Concentration	Batch dependent - please inquire should you have specific requirements.
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.