



Bovine IgG [HRP] (DAGA-3141)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|-------------------------|---|
| Product Overview | Bovine Immunoglobulin G (HRP) |
| Purity | Bovine IgG (HRP) was purified by delipidation, salt fractionation and ion exchange chromatography followed by dialysis. |
| Conjugate | HRP |
| Applications | User optimized |
| Format | Lyophilized from 0.01M Na ₃ PO ₄ , pH 7.4, with 0.14M NaCl, and 10 mg/ml BSA. Immunoglobulin and Protease free. |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Storage | Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles. |